

Where water becomes the focus

- ✓ Adapted to safety standards.
- ✓ Sink workstations with water taps and fittings, draining racks, soap and towel dispensers and eyewash stations for first aid.



Top-mounted basin

- All-around marine edge; larger than the worktop opening

Built-in basin


- Set into the worktop

Welded-in basin

- Welded seamlessly into the worktop

Underslung basin

- Fastened underneath the worktop





	Materials	Widths	Depths	Configuration options
Basins	Ceramic, stainless steel, polypropylene, epoxy	250 – 895	100 – 500	Top mounted, flush mounted, underslung
Sinks	Ceramic, stainless steel, polypropylene	1200	600 750	Basin left/centre/right, with/without drainboard
Front sinks	Ceramic, stainless steel, polypropylene, epoxy, Trespa® TopLab® Plus	1350	680	Basin central, drainboard left and right
		1500	680	
		1650	680	
		1800	680	
		1950	680	

Dimensions [mm]



Basin material / worktop material	Ceramic	Polypropylene	Stainless steel	Epoxy resin
Solid ceramic	+	-	-	-
Composite ceramic	+	-	-	-
Composite glass	+	-	-	+
Epoxy resin	-	+	+	+
Trespa/HPL	-	+	+	+
Melamine	-	+	+	+
Polypropylene	-	+	-	-
Stainless steel	-	-	+	-

Dimensions [mm]