

Service cells to meet diverse needs

Floor-mounted: Positioned behind the worktop or separate supply


Bench-mounted: Fixed laboratory bench and cell unit for a greater work area



For the best possible space utilisation

- ✓ Attractive appearance and translucence
- ✓ Workplace structuring from vertical arrangement of services
- ✓ Cell cabinets and shelves provide additional room for equipment and utensils





	Width				Height above top edge of floor			Height above worktop			Depth
Service cells for sitting workstations	900	1200	1500	1800	1320	1620	1920	570	870	1170	150
Service cells for standing workstations	900	1200	1500	1800	1320	1620	1920	420	720	1020	150

Dimensions [mm]