


Utensils always ready to hand

- ✓ Customized storage space as a mobile cabinet, frame base cabinet or solid with a steel plinth
- ✓ Flexibly and individually configurable
- ✓ C-frame or A-frame
- ✓ Worktops matched to your requirements



	Width				Height		Carcass depth	
With steel plinth	450	600	900	1200	720	870	350	500
With steel plinth for front sink	300	-	-	-	720	870	-	500
With castors	450	600	900	1200	580	730	-	500
For installation in frames	450	600	900	1200	480	630	-	500



Dimensions [mm]