

# CERTIFICATE

It is hereby confirmed that the

**Fume Cupboard  
For bench type 1500 mm**

of the manufacturer

**Köttermann GmbH & Co KG  
Industriestraße 2 - 10  
31311 Uetze/Hänigsen  
Germany**

was tested according to

**EN 14175  
part 3**

at the total volume flow rate of **490 m<sup>3</sup>/h**,  
with a nominal face velocity of **0.2 m/s**  
With this exhaust volume flow rate the requirements of  
- **BG Chemie, Germany, issued 29.07.2003**  
- **XP X 15-206, France, from Jan.2005**  
are completely fulfilled

The details and data of the complete test procedure  
are content of the type test report

**FCS/EN/184/2013**

which is given to the manufacturer on  
24. June 2013

Tested by:  
Tintschl BioEnergie und  
Strömungstechnik AG  
Goerdelerstr. 21  
91058 Erlangen

*Sulay Tröger*  
.....  
Dipl.-Chemie. Sulay Tröger  
Erlangen, 26.06.2013

